



Please read SAFETY INSTRUCTIONS carefully before using the chair.





Titan Chair Introduction to the Titan Jupiter

- 1. The intelligent 3D technology will deliver a consistent soothing massage from your neck to down to your hamstrings. Adopting traditional Oriental massage methods, you will come to experience a soothing kneading, pressing, tapping and shiatsu style massage that will invigorate the senses and relax the body.
- 2. With acupuncture point tracking innovation. It will deliver a massage that will massage the key areas while adjusting to different body type.
- 3. To enhance the massage experience the backrest is fitted with a warm soothing heating pad
- 4. The Jupiter massage chair is equipped with a highly advanced L track system that reaches neck to the hamstrings.
- 5. One simple touch of the remote, recline back into a comfortable zero gravity position that promotes an even distribution of weight throughout the back.
- 6. The forward sliding space saving technology is one of the great features that greatly eliminate the awkward space required to recline. The backrest and seat frame glides forward as you recline the chair.
- 7. With advanced wireless Bluetooth technology, mobile or tablet device can be synced with the massage chair, to answer calls or listen to your favorite tunes.
- The foot massager has a rounded points that scrape the bottoms of the feet in 8. combination with air bags that apply a squeezing/kneading action.
- 9. In the calf region the airbags will deliver a pulsating action that that will squeeze the and release, refreshing the legs.
- 10. For the arms the same type of airbag massage is delivered. Simply slide your arms into the slots soon as you select the massage of your preference.
- 11. Multiple airbags are thoughtfully placed throughout the chair to provide a fully body massage experience.
- 12. The legrest and foot parts can be adjusted to any angle that you needed, It can extended and retracted within the range of 15cm, to accommodate varies heights.
- 13. The high resolution LCD screen allows the user to see what the chair is performing and where the massage roller heads are located.

- 1. Thank you for selecting the highly advanced Jupiter massage chair by Tltan. 2.Please read this manual carefully before using the massage chair. 3.Please keep this manual in a safe place for further reference. 4.Please pay attention to SAFETY INSTRUCTIONS.
- \checkmark Note: We reserve the right for design modification. It is subject to change without notice. The pictures in this manual are only for reference. Please refer to the actual product.

Preparatior

Massage Operatior

Content	
Safety instructions	
External structure parts	
Internal structure parts	
Mechanical components diagram	
Mechanical components list	
Installation	
The adjustment before massage	
Manual control	
Operation instruction & function illustration	
After massage	
Clean and maintenance	
Malfunction guide	23
Technical dat a	24

Please keep this manual in a safe place for future reference.



Safety instructions

- Please read this manual carefully to make sure the installation and operation of the massage chair is correct.
- Safety Instructions are provided in in this manual with Warning and Caution to ensure the safety of the user.

Warning and Caution are highly important so please be aware of all safety warnings.



Symbol meaning

Preparation	
	l

This sign indicates forbidden action.

This sign indicates self disassembly is forbidden.

operation of cable.

This sign indicates the correct

This sign indicates you must follow instruction to operate.

The following should be kept for any further references.

Warning person with Osteoporosis are not permitted to us • People with heart disease including implanted other medical electronic device are not permittee For pregnant women or people with any medical please consult doctor before using. Patients with a high fever are forbidden to use the Perceptual disabilities patients are forbidden to People who are elderly should consult with their Children must consult with their doctors prior to Do not sleep in the chair. THis may lead to serious Don't use the chair when your body is wet to pr Make sure there are no kids or pets or other object (under, behind or in front of)when the chair is in ope • During operation, if you feel any abnormal discomf the operation and immediately consult your physi Please do not use the chair to massage above you • This chair is intended to be be used in the sit down • Do not use this chair while intoxicated. If you have any neck or back issued you must get ap doctor. This chair produces a minor electrical field, therefore metal plating or medical devices should consult y Insert the plug of the chair into correct power out plug is completely inserted into the socket to circuit or fire. • Unplug the power when not in use, stop all mass and turn off before unplugging. • It is prohibited to use the massage chair with an coverings which may lead to damage or injury. Massage cushions of the chair are not permitted to removed to prevent injury. The machine is forbidden for use if no one is sitt \mathcal{O} • Be sure the power cord or any wires are in good c Don't use other electric power except rated voltage electric shock or any breakdown. Any broken, re-process, knotted, twisted power of under chair is not allowed for use. If the chair has been exposed or made contact with any float and unplug the chair from the outlet until it has been inspec by a certified technician. • If power cord or plug is damaged, please don't replace contact us or appointed dealer for replacement. To avoid electric shock or bodily injury please don't op (\mathbb{R}) any part of the chair. The chair should not be handled unless approved by the manufacturer. • Check the upholstery or cloth material on the back a damage or extreme wear. If broken, even though it is using the chair, unplug the power and have it repaired major damge.

Safety

y instru	uctions	2
se this product. pacemaker or ed. al condition, this machine. use this device. doctors prior. using. us injury. revent injury.	People who are not suitable to use the chair	
jects around eration. Ifort, please stop sician. Dur neck. postion. approval by your fore anyone with your doctor.	Accident/injury	Preparation
tlet. Make sure the prevent short sage functions	Fire and electric shock	
ny blankets or to be folded or	Accident	
tting in the chair. condition. age to avoid fire, cord or pressed	Fire and electric shock	
fluids STOP using ected and repaired peen or disassemble d by the customer and other areas for small, please stop d to prevent from		

Safety instructions

3

Preparation

	Caution	
	• The chair must be placed on a flat even surface.	To avoid sudden falling.
	 The backrest of the chair when upright must be placed a minimum of 10 cm from the wall. 	To make sure the back lying.
	 The maximum suggested operating time is 20mins. A maximum of 5 minutes for the upper, mid and lower back is advised to avoid injury. 	Too long time massage is bad for your body.
•	 Before sitting on the chair to avoid injury please inspect everything to confirm that the massage heads are in the correct position. Remote control should be put in the pocket or in the bracket. Do not put the remote control anywhere besides the suggested locations to avoid damage. 	To avoid accidental injury.
U	• Please make sure that the power switch is at the "OFF" position before inserting the plug into the socket.	
	• Turn off all remote before removing plug from outlet.	To avoid the accident when turn on the power.
	• Make sure the plug is completely insert into outlet.	To reduce short circuits and any fire hazard.
	 Checking for any damage on wire and make sure there is no fluids are near or on wire. 	
	 Do not tamper or remove the ground wire The plug must be inserted into socket with the ground prong in tack and in good condtion. 	To avoid electric shock in any leakage or breakdown.
\bigcirc	 Don't press down on the massage heads with hands or step on the massage heads with feet. this will damage the chair. Don't put hands between the working massage heads. Don't put hands or feet in the gap of mechanical parts. Don't drop or embed anything into the cracks of massage chair. Don't stand or sit or put anything on the backrest, armrest or leg-rest. Don't stand on the chair during operation. Don't put lout the plug or turn off the power during operation. Don't put your hands or head between leg-rest or beneath it. Don't use the massage chair while wearing hair accessories. Remove everything in your pockets or hard objects that may cause injury. Don't apply additional force on the backrest when sitting in the chair. When using the this product, you must be properly clothed to prevent injury or irritation to the skin. 	To avoid accidents or injury
	 Keep the chair away from heat, fluids, sharp things and corrosive, inflammable, explosive products. To use electric blanket or other heating products on the chair is forbidden. Don't place near any corrosive or or flammible products Don't put the chair in a small space places or use in areas with poor air circulation. Keep the chair away from stove or other heating products. Do not expose the chair to sunlight or areas with temperatures above 85 degrees Fehreinhiet. Don't use the chair outdoors or sunrooms. 	To avoid fade, discolor on upholstery or hardened, deteriorated on PU.

	Do not place this chair in areas that have above normal humidity	
S	 Bo not place this chain if aleas that have above normal numbers levels. Keep away from fluids and must be dry at all times. The power cord must be the original one provide. If not contact the manufacturer. Do not use until approved by manufacturer or certified electrician. 	To avoid leakage or electric shock.
ļ	• It is imparitive that all electrial wires are in good condition and that there are no visiable cracks or exposures of the wires. Be sure the power cords are all firmly inserted into the chair and the wall outlet. When using the chair, no part of your body should be wet, especially the hands.	To avoid electric shock.
	 Unplug the powercord form the outlet before cleaning chair. When the chair is not being used unplug from the outlet. 	-
	If the chair has any burning smell immediately unplug the power.	To avoid electric shock or fire harzard.
Z	• When removing the plug from the outlet, firmly hold form the head of the plug. Do not pull from the wire portion.	
	• Please turn off the power switch and unplug the power after massage.	To avoid accidents
	Unplug the power immediately when electricity goes out	To avoid damage or accidents when power on again.
	ding instruction * The chair is equipped with a cord having an grounding conduct plugged into an appropriate socket that is properly installed and matched with the plug.	tor. The plug must be
Repla	ce fuse	
\Rightarrow	When the fuse is burned out, please turn off the switch and unplug the power. If the chair is not powering up, please check the fuse by removing plug cover shown on the right Using a same size fuse with the same rating to replace the old one, then cover the lid. Or you can	
	contact the manufacturer for a certified electrian.	

Safety instructions 4





External structure parts

5





Mechanical components diagram

Explosion Diagram And the second second \$

Preparation

		chanical c	omponents	STISC	
NO.	Name	Q' ty	NO.	Name	Q' ty
1	Remote control	1	22	Shoulder pads	1
2	Right side fender	1	23	Pillow	1
3	Air pump assembly	1	24	Right armrest assembly	1
4	connecting rod	1	25	Seat frame surface assembly	1
5	Legrest & footrest	1	26	Seat support bar	1
6	Front fender	1	27	Castor	2
7	Left side fender	1	28	Power case assembly	2
8	Electric magnetic valve	1			1
9	Backrest electric actuator	1			
10	Bottom seat frame assembly	1			
11	Footrest electric actuator	1			
12	Sliding rail	2			
13	Top seat frame assembly	1			
14	Left amrest assembly	1			
15	Electric magnetic valve	2			
16	Shoulder massage device	1			
17	Massage machine	1			
18	Back cover assembly	1			
19	Back frame assemly	1			
20	Backrest plastic front cover and seat plate assembly	1			
21	Back and seat cushion	1			

Mechanical components list

		,	
	2		
/			

Preparation	

8





Step 1. How to move the chair

①Cut the straps and unpack the box to access the base of the chair. Remove all pads and cushions from the base and lift the chair from the points indicated in image below.
With 2 or more people, lift carefully and move the chair to the desired location.
②Be sure to keep all the cushions, remote and power cord together so that nothing is lost.

(as the picture shows)



the armrest carton box. Remove the armrest slowly and carefully. We recommend that you lay the armrest face up on a blanket or carpeted floor to avoid scratches and tears. When holding the armrest, handle by grabbing on the indicated areas (image on the right).



(4) Take out the legrest from the legrest box. remove from the box by grabbing the indicated locations.



Installation (2)

Step 2. Massage chair position return

1.Once the chair base and backrest are removed from the box, insert the power cord into the chair and then into the outlet 2. Turn on the power switch on the back of the chair and the backrest will raise up as the seat base levels flat. 3.Once the chair is in the upright position turn off the chair and proceed to next step of the install. (Refer to step 7, after the assembly)



Step 3. How to assemble the legrest

Install the legrest to the main chair body. 1. Unscrew the bolt from the hanger of legrest . 2. Please insert three air tubes of the legs into the air taps of the main body. Connect the air tubes according to the color code. 3. Above the air connections, insert the power connector into the socket firmly. 4.Lift up the legrest and align the axis into the hanger of main body, then put down the axis cover. 5. Screw on the bolt tightly.







Massage chair position return state

Installation





The adjustment before massage

Environment Check environment Space size for using the chair ∕!∖ Caution • At least 10cm from wall. Don't use the chair in wet places like near swimming pool, bath room to avoid leakage or electric shock. Don't use under sunlight or near stove a star or other heating place to avoid upho-Istery deteriorated. Please use chair in flat place to avoid chair falling or noise or other unexpected trouble. • Check places in front Check the areas when back and leg extending, • At least 30cm in front of the leg. make sure no people, pets or other things Check surroundings around. When the backrest or legrest reclines or extends, please pay attention to: (1) Check any person or things in every direction. NVD $\nabla \nabla P$ (2) Check the space size. Ø Check power cable, plug (1) Clean the dust on plug on time. (using dry cloth). (2) No damage, re-do, bend, stretch, tie or press with heavy thing(especially can't be pressed by the chair) on power cable. Check

ow to move massage chair
Don't move chair when it's working in case of accidents or damage
How to move Bind the cables and put them in the seat to avoid damage when moving. If you plump the chair, it may damage the inner parts of the chair . Don't use the castor if the floor is wooden or other easy damaged material. It will need 2 or more persons to lift the chair. Make sure the power is off and the plug doesn't touch the floor(It's better to let the chair in standing status to avoid theback crashing the floor)
Lift up the legrest to make the gravity fall on the castor. Then push chair to destination and set it down.
Please hold the footrest to move the chair as shown in the picture in case any scratches on the leatherette

Massage operation

The adjustment before massage



e adjustment

e make sure the age head in the art before sitting e chair.

id injure caused by position of the ge head.



Incture point tracking detecting

doing acupuncture point position ting, if massage head can't detect the ve shoulder height position, the actual der height position was default to the e part of back.

position: Your back the chair back, the lings to the pillow.

mote control

'beep" sound after

ncture point tracking detecting ed,when micro-adjustment is needed, ress up/down key on the remote

for the massage head detecting your nctuer point postion not right,otherwise on't get good massage.

suggestion



 When you feel the intensity isn't strong enough,take off the pillow or cushion.

e the chair is heavy, it may ruin or in long time setting.

put a blanket or oft things to your floor. Massage operation



Manual control



Manual control

Massage operation



interface. There are 10 massage programs to select from: power, sleeping, vitality mode, relax, comfortable mode, pain relief, muscle strengthen, blood circulation, focus on waist, stretch.

All pre-set massage programs will conduct a computer body scan. Soon as the computer body scan is complete, choose any one of the 10 programs, shoulder height detection is not needed any more.

Back/seat angle can be adjusted automatically: During auto massage, the back/seat angle will be automatically changed.

Manual control 18 LCD display graphic description Manual Massage: If the computer Manual body scan doesn't complete, Manual please restart the process again. In 3D Feel Manual Neck&shoulder manual mode the computer body Position Position Kneading Shoulder&back scan will activate unless "point" Fixed Point Fixed Point Back&waist Tapping massage is selected. Knocking Waist&sea Shiatsu Overall 1.7 massage methods: shiatsu, Mix 1 kneading, tapping, knocking, shiatsu Mix 2 (rolling), mix 1(Knead+Tap) and mix 2(Knead+Knock). 2.5 positions for choice: Select by the navigation keys Select by the navigation keys Neck and Shoulder Mode: massage neck and shoulder Shoulder and Back Mode: massage shoulder and back Back and Waist Mode: massage back and waist Waist and Seat Mode: massage waist and seat Whole Body Mode: massage the whole body

During 3D Feel mode, only three positions are available, neck and shoulder mode, shoulder and back mode, back and waist mode.

- 3. Fixed Point
- . The preset massage program default mode is whole body.
- 2. Select massage method before choosing the position or fixed point.
- 3. Select Manual mode, all air pressure massage will turn off.





Air Massage: Ten air massage modes for you to choose, Shoulder+Arm, Waist+Pelvis, Leg+Foot, Head, Shoulder+ Foot, Arm+Pelvis, Waist+Foot, Upper Body, Lower Body, Whole Body Mode. Once air massage is activated, the mechanical massage of the back massage machine shall stop.

Operation instruction and function illustration



Massage operation



After massage

Clean and Maintenance parts into neutral detergent, it and clean. (FIGO) h water-dipped but d cloth. dry. ne controller th soft cloth after squeezing. dry. nt, diluent or re prohibited. **I** Controller About cleaning leaning, unplug the electricity, uch plug with wet hands. re hands are dry when plugging cket. to follow the instruction you are otential serious injury.



Synthetic leather	Plastic
General cleaning Clean with a clean soft and dry cloth.	 Dip cloth squeeze Wipe with squeezed Allow to c
Forbidden interest poster If upholstery is dirty	Clean th
Dip soft cloth into 3-5% neutral detergent solution, squeeze it then clean the stain on the surface gently.	Wipe wit
 Wash cloth, squeeze it and wipe with detergent. 	Allow to e Propellan alcohol a
Iclean it softly with dry cloth.Allow to dry (air blower is prohibited)	Forbidden
Cloth material	
 Dip cloth into weak neutral detergent, squeeze it and clean the surface & fabric. (Diluent, gas, alcohol are prohibited) Brush the dirty parts with weak neutral detergent. (avoid excessivebrushing) 	 Before cl don' t tou Make sur to the source

23 Malfunction guide

Trouble	Possible reason	Troubleshooting
	Check if the control is on or not.	Open control switch
	A massage function was not selected	Choose massage function
The measure sheir door	check if the power cord is firmly inserted into the outlet and the chair.	Make sure plug and socket are connected
The massage chair does not work after being powered on.	Wire or plug has been damaged.	Please ask manufacturer or appointe or authorized technician to inspect an repair.
	The fuse burned out.	Replace fuse of same specification.
	Something wrong with inner circuit	Please ask manufacturer or appointed certified technician to inspect and repair.
There is beep sound from the inside while working.	Sound is originating from air pump or motor.	This sound is normal once the chair has broken in.
Different height of the left roller and right roller.	The chair is working fine but slightly off level.	The arm or spring is bent or broken
Sound suddenly becomes	Check if the chair has been running for more than 60minutes	Turn off the power to have the massager rest for over half a hour before next usuage.
much louder while working.		
Massage function suddenly stops while working.	Suddenly power-down timer has e×pired.	Turn off the power to have th massager rest for over half a hour before next turn working.
Neither the backrest nor the legrest can be lay down.	Check if there are any obstacles; the chair overloads; for the sake of safety, sometimes the massage machine may	Remove the obstacles; turn off the power to have the massage rest for over half an hour before using again.
The chair can't return to its original position.	automatically stop working. Excessive operation; the timer has	Please ask manufacturer or appointe or authorized maintenance man for checking and replacement.
Electric wire or plug is	e×pired The chair may be overused	Turn off the power to have the massager rest for over half a hour before using again.
e×tremely hot.	There may be an electrical issue	Stop using and call your manufacturer immediately. Be sure to unplug power

Technical data

Name	Massage Chair	Power
Model	RK7909	NW/GW
Input voltage	□ AC120V, 60Hz □ AC220V, 50/60Hz □ AC240V, 50Hz	
Safety design	Class I	Packing size

/!\

Note: If the troubles are not caused by above reasons, please turn off the power switch and cut off the power, then contact the manufacturer or service agent to repair your chair.

Technical data 24 200 W Main body :78/102 kg Armrest : 15.5/22kg Legrest & footrest : 15/17.5kg Main body:1510×790×1030 mm Armrest : 1245×410×655 mm Legrest &footrest :530 × 445 × 540 mm $(L \times W \times H)$ Others